

# User Manual



Conpreva<sub>2</sub>



# Content

<b>1 IMPORTANT</b> .....	<b>4</b>
Intended use .....	4
Medical application areas .....	4
<b>2 SAFETY AND WARNING INSTRUCTIONS</b> .....	<b>5</b>
Explanation of graphic symbols and product markings .....	5
Safety instructions .....	7
<b>3 SCOPE OF DELIVERY</b> .....	<b>13</b>
Optional accessories .....	13
<b>4 TECHNICAL DESCRIPTION</b> .....	<b>15</b>
Device layout & accessories .....	15
Overview of functions and signals .....	16
Environmental conditions during use .....	17
<b>5 TRAINING PREPARATIONS</b> .....	<b>18</b>
Preparing Electrodes .....	19
Preparing Sensors .....	19
Preparing Pressure Measurement .....	19
<b>6 LAUNCH EXERCISE</b> .....	<b>20</b>
Therapist & Patient Selection .....	20
Program selection .....	20
EMG-Setting .....	22
Electrical current .....	23
Pressure Setting and Zero Adjustment .....	24
Exercise process .....	24
EMG-biofeedback training .....	25
Electrical stimulation training .....	26
Combined training (Biofeedback + electrical stimulation) .....	27
EMG-triggered electrical stimulation .....	29
Pressure measurement .....	31
Pressure measurement + 1-channel EMG .....	31
Pressure measurement + 2-channel EMG .....	31
Pressure training .....	31
Finishing a therapy session .....	32
<b>7 MAINTENANCE</b> .....	<b>33</b>
Maintaining sensors .....	33
Maintaining electrodes .....	33

Maintaining the device.....	33
Replacing batteries.....	34
Safety checks .....	35
Troubleshooting.....	35
Storage and transport.....	37
<b>8 ENVIRONMENTAL PROTECTION &amp; DISPOSAL ...</b>	<b>38</b>
<b>9 TECHNICAL INFORMATION.....</b>	<b>39</b>
Name of the medical device .....	39
Model .....	39
Display .....	39
Connection ports .....	39
Modes of operation.....	40
EMG-Range.....	40
Output.....	40
Pulse pattern .....	41
Pulse group pattern (medium frequency) .....	41
Exercise phases .....	41
Power supply .....	41
Safety features .....	42
Operating conditions .....	43
Dimensions and weight .....	43
Storage of exercise data .....	43
Application recommendation .....	44
ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES	45
<b>10 WARRANTY .....</b>	<b>49</b>
<b>11 ADVANCED INFORMATION FOR PHYSICIANS AND THERAPIST .....</b>	<b>50</b>
Therapy Settings .....	52
System Settings.....	54
Administration.....	58

5th edition, October 2017

This instruction manual contains copyright protected information of Haynl-Elektronik GmbH and third parties. Each rendition or utilization outside of the borders permitted by copyright is inadmissible without previous written permission on the part of Haynl-Elektronik GmbH.

Haynl-Elektronik GmbH reserves the right to make improvements and changes without prior notice. All data and information is reviewed periodically and any changes applied will be implemented in the next available version of the instruction manual.

The Bluetooth® brand name and logo are property of Bluetooth SIG, Inc. All usage of this trademark by Haynl-Elektronik GmbH is covered under a license agreement. Other trademarks may contain property of other owners.

# 1 IMPORTANT

PLEASE READ ALL INFORMATION PROVIDED IN THIS INSTRUCTION MANUAL BEFORE USING THE CONPREVA<sub>2</sub> MEDICAL DEVICE OR ANY OF ITS ACCESSORIES.

THIS INSTRUCTION MANUAL IS AN INTEGRAL PART OF THE MEDICAL DEVICE.

PLEASE STORE THIS MANUAL WITH THE MEDICAL DEVICE.

Do not use Conpreva<sub>2</sub> if you do not fully understand any aspect of operating the device. Please refer to a qualified supplier for assistance.

## ***Intended use***

The electrotherapy device Conpreva<sub>2</sub> is a class IIa medical device. It is intended for the electrotherapeutic training of pelvic floor musculature.

## ***Medical application areas***

The Conpreva<sub>2</sub> therapy device is designed for the treatment of patients afflicted with stress incontinence, urge incontinence, fecal incontinence or mixed incontinence, bladder sphincter dyssynergia and for postoperative training following prostatectomy or imperforate anus.

## 2 Safety and warning instructions

### *Explanation of graphic symbols and product markings*



Indicates a major hazard that poses a high risk of death or serious injury if not avoided.



Indicates a significant hazard that poses a medium risk of death or serious injury if not avoided.



Indicates a minor hazard that poses a risk of medium to minor injuries if not avoided.



Consult instruction manual



Manufacturer



Serial number



Application part Type BF



„On“/„Off“

CE 0123

CE-marking and identification number of notified body involved



Caution – follow instructions



Protect medical device from moisture



Protect device from direct sunlight and heat



Do not dispose of medical device along with normal/household waste



Not suitable for pacemaker wearers



Prevention of danger/risk



Article number



Serial number



# Safety instructions

## DANGER



### **Risk of interference with cardiac pacemakers**

- Conpreva<sub>2</sub> should not be used by patients fitted with a pacemaker unless with explicit permission given by a medical specialist.

### **Danger to patients with cardiac damage**

- Conpreva<sub>2</sub> should not be used by patients with cardiac damage unless with explicit permission given by a medical specialist.

### **Danger to children and untrained users**

- Store device at a secure location.

### **Risk of damage to implanted medical products**

- Conpreva<sub>2</sub> should not be used by patients with medical implants unless with explicit permission given by a medical specialist.

### **Danger to patients during pregnancy**

- Conpreva<sub>2</sub> should not be used by pregnant patients unless with explicit permission given by a medical specialist.

### **Deviation of electrical output or electric shock caused by liquids inside the device**

- Do not expose the device to high humidity. Do not clean the device with water or other liquids. Do not use the device outdoors.

### **Deviation of electrical output by unauthorized repair attempts**

- Do neither open nor attempt to repair the device.

### **Risk of igniting highly flammable substances**

- ➔ Keep the device away from easily flammable materials.

### **Usage by multiple patients:**

#### **Health risks due to contamination, chemical residues and changes of functional performance**

- ➔ Requires validated processes in accordance with latest legal requirements to prepare the device for use by a different patient.

#### **Risk to health (heart, head, genitalia, spine, neck) due to electrical stimulation**

- ➔ Only use the device as per the instructions at all times.

#### **Damage to health caused by incorrect use of device**

- ➔ Only use the device as per the instructions given in this manual at all times. Always store this manual together with the device.



## **WARNING**

#### **Risk of skin burns due to high-frequency surgery equipment**

- ➔ Simultaneous use of a surgical device employing high frequencies and this device may cause burns at contact points between skin and electrodes.

#### **Risk of infection**

- ➔ Any accessories may only be used by a single patient.

#### **Risk of flammable/explosive substances being ignited**

- ➔ Do not operate the device near flammable/explosive substances or flames of any type.

**Deviation of electrical output due to ambient conditions**

- ➔ Ensure that the recommended environmental conditions for using and storage of the device are not exceeded.

**Risk of infection for patients designated to use sterile products only**

- ➔ Do not use the device if sterile environmental conditions are required for a patient.

**Danger due to erratic or involuntary movements**

Do not use the device with patients suffering from seizure disorders (e.g. epilepsy), increased epileptogenicity or severe cognitive impairments unless with explicit permission given by a medical specialist.

**Risk of electrical shock or injury caused by non-approved accessories**

- ➔ Only use approved accessories with the device.

**Risk due to improper use**

- ➔ Do not connect Conpreva<sub>2</sub> with any accessories or medical devices other than those approved by the manufacturer.

**Risk of electrical shock from electricity grid**

- ➔ Do not connect the device, any of its parts or accessories with the electricity grid.



## **CAUTION**

### **Deviation of electrical output in the proximity of therapeutic shortwave or microwave devices**

- ➔ Do not operate the device near therapeutic shortwave or microwave devices.

### **Pain or physical discomfort during treatment**

- ➔ Stop using the device immediately, consult the medical specialist in charge of the treatment and notify the manufacturer.

### **Danger of injury due to damaged accessories**

- ➔ Avoid pulling on connection cables or applying excessive force to any part of the accessories. Do not use any damaged accessories.

### **Danger of injury due to damaged device casing**

- ➔ Avoid damaging the device casing by applying excessive force. Do not drop the device. Do not use damaged devices.

### **Danger of injury by electrical shock at the start and end of a therapy session**

- ➔ Only turn the device on after all accessories have been connected. Do not remove any accessories before the device has been switched off.

### **Danger of injury by electrical shock when changing batteries**

- ➔ Change batteries only after device has been switched off and all accessories removed.

### **Danger of injury by electrical shock**

- ➔ Do not switch device on while battery compartment is open.

**Danger of tissue damage around implants**

- Do not place electrodes near implants containing conductive materials.

**Danger of injury due to surgical wounds, skin injuries or skin irritations**

- Do not place electrodes on or near any skin injuries, skin irritations, open wounds or fresh surgical wounds.

**Danger to health due to infectious diseases or fever**

- Do not use device.

**Danger of injury while operating machinery**

- Do not use device.

**Danger of injury while driving a vehicle**

- Do not use device.

**Danger of injury while under the influence of alcohol or medication which may impair the ability to operate machinery**

- Do not use device.

**Danger of injury by high electrical currents in case of sensitivity disorders**

- Configure the intensity level of the electrical stimulation together with the medical specialist in charge of the treatment.

**Danger of health impairments for patients suffering from polyneuropathy or peripheral nerve lesion in the region of the pelvic floor**

- Do not use device.

**Danger of triggering allergies**

- Do not use device or its accessories if patient is allergic to gold or plastics.

**Danger of bleedings due to haemorrhoid related illnesses**

- ➔ Do not use device or accessories unless with explicit permission given by a medical specialist.


### 3 Scope of delivery

 <b>WARNING</b>
<b>Danger of infection</b>
→ Any accessories may only be used by a single patient.

The Conpreva<sub>2</sub> therapy device is delivered as a set containing the following items:

- 1x Conpreva<sub>2</sub> (HPS12)
- 4x AA - Battery
- 1x Instruction Manual

#### *Optional accessories*

 <b>WARNING</b>
<b>Danger from third-party accessories</b>
→ Only authorized third-party accessories (medical products) with CE-mark may be used for the intended purpose within the limits of the medical application;
→ Only use third-party accessories mutual compatibility between the products with respect to manufacturer guidelines, the required steps during use, manufacturer guidelines for use and packaging has been reviewed and a statement of such is made available; The actual use is to be overseen and verified in a suitable way.

A number of optional items for use together with the Conpreva<sub>2</sub> are available. Some of them are required in order to perform certain treatment exercises while others enhance the general function of the device.

Any optional accessories must fulfill all relevant requirements listed in this manual and the regulatory requirements demanded of medical products.

**Electrodes may only be connected via K1 (MED) lead cables and K3 (MED) lead cables with a device.**

The following items are available:

- Rectal-sensor (to connect with K1 (MED) cable)
- Vaginal-sensor (to connect with K1 (MED) cable)
- Type I: 4 sticking electrodes 4x4cm angular
- Type II: 24 sticking electrodes  $\varnothing$  2cm round
- Sticking electrodes of different shape, size or placement as assigned by medical prescription to be connected with K1 (MED) or K3 (MED) cable
- Rectal pressure sensor with Luer-Lock-adapter (female)
- Vaginal pressure sensor with Luer-Lock-adapter (female)
- Three pole cable K1 (MED), 2mm plug
- Three pole cable K3 (MED), 1.5mm plug
- Contact gel
- HPS-Software for creating new exercises and review of exercise data
- Medical USB-cable for PC connection
- 4x NiMH rechargeable batteries (2500 mAh)
- 1x external battery recharger
- Power supply unit (6V  $\overline{\text{---}}$  / 1A; 1.3mm plug;  $\varnothing$  3.7mm max) conform with **EN60601-1:2006** and **EN60601-1-2:2007**



## WARNING


### Danger due to high current density

- ➔ When using electrodes with a surface area smaller than 5 cm<sup>2</sup> (square of 2.2 cm or circle of 2.5 cm) device operators should be very attentive during device usage



# 4 Technical description

## Device layout & accessories




# Overview of functions and signals


## Acoustic signals

- )) Single short beep: button pressed
- )) )) Two short beeps: batteries low, device will turn off


## Optical signals and icons



- Muscle activity in  $\mu\text{V}$
- Icon for „Continue“ (Button A)
- Icon for „Submenu“ (Button B)
- Bar graph showing muscle activity
- Icon for „Increase range“ (Button C)
- Icon for „Decrease range“ (Button D)
- Current range in  $\mu\text{V}$



- Electric current intensity in mA
- Icon for „Continue“ (Button A)
- Bar graph showing intensity of electric current
- Icon for „Increase current“ (Button C); appears only after Button A has been pressed
- Icon for „Decrease current“ (Button D)



- Caution! Current interrupted; check cable and sockets
- Icon for „Increase current“ (Button C)
- Icon for „Decrease current“ (Button D)

## ***Environmental conditions during use***



### **CAUTION**



#### **Danger to operator/user due to defective device**

→ Medical device may only be operated at an ambient temperature between 5°C and 35°C.



→ Protect device from exposure to moisture.



→ Protect device from strong sunlight and strong sources of heat.

→ Do not use device in dusty environment.



→ Medical device may only be operated at a relative humidity of 30-75%.



→ Medical device may only be operated at a relative pressure of 700-1060 hPa.

→ Therapy device may only be used in extreme temperature fluctuations (for example in winter) after sufficient acclimatization (approx. 1h)

## 5 Training Preparations

Please choose a quiet and relaxed environment for therapy sessions. You should be able to sit comfortable or lie down. Avoid the proximity of TV sets, loudspeakers, fluorescent lamps or similar devices. The Conpreva<sub>2</sub> medical device is equipped highly sensitive sensors. As such strong electromagnetic fields may cause significant errors during measurements despite appropriate shielding.

Before starting a training session, the device needs to be turned off.

The intensity levels of electrical stimulation should be chosen in accordance with the medical specialist in charge of the treatment.

Please read the passage *Preparing Electrodes* or *Preparing Sensors* depending on whether electrodes or a sensor have been recommended for the treatment.

The treatment should never feel unpleasant or painful; however, a clear electrical tingling sensation should be felt.

## ***Preparing Electrodes***

- Thoroughly clean the skin areas where the electrodes are to be placed. The skin needs to be dry and hairless. Use contact gel to improve conductivity and adhesion between skin and electrodes.
- Only apply electrodes to skin areas as recommended by the medical specialist responsible for the treatment.
- Connect the device plug of the connection cable with the *CH1 and/or CH2 socket* on the device.
- Connect the electrode plugs of the connection cable with the electrodes.
- Turn the device on by pressing the *On/Off Button*.

## ***Preparing Sensors***


- Insert the sensor rectal (men/women) or vaginal (women).
- Connect the plug of the connection cable with the *CH1 and/or CH2 socket* on the device.
- Turn the device on by pressing the *On/Off Button*.

## ***Preparing Pressure Measurement***

- Connect the plug of the pressure sensor with the *CH2 socket* on the device.
- If a combined measurement (pressure + EMG) is to be conducted, please follow the instructions under *Preparing Electrodes* or *Preparing Sensors*
- Turn the device on by pressing the *On/Off Button*.

## 6 Launch Exercise

### *Therapist & Patient Selection*



After turning the device on by pressing the *On/Off Button* it will show a boot screen. Once the device has finished booting up its screen will display a selection for therapists and patients.

The name of the currently selected therapist is displayed in a blue box at the top of the screen. The assigned patients of the current therapist are shown in numbered green boxes taking up most of the screen.

The buttons  and  next to a therapist's name can be used to change the current therapist. To select a patient simply press the green field with their name.

There can be up to six different therapists. Each therapist can be assigned up to 6 patients. Each patient in turn can have up to 6 different exercise programs available

### *Program selection*

After the patient selection the screen will show a selection of programs available for this patient. The program selection may cover two screens if four or more programs are available.

Press *Buttons A to C* to start the treatment program displayed next to it. Press *Button D* to switch between

screens of the program selection. Only one program may be started at a time.

The background color of a treatment program indicates its mode of operation:


- Blue: EMG-Biofeedback



- Red: Electrical stimulation




- Green: Combined training (EMG + stimulation)



- Violet: EMG-triggered electrical stimulation




- Cyan: Pressure measurement




- Cyan: Pressure training



- Cyan: Pressure + 1-channel EMG




- Cyan: Pressure + 2-channel EMG



Details concerning each program are displayed below their names. This includes the length of different parts of each program such as Tense (AZ), Relax (EZ), Stimulation (SZ) and Break (PZ) phases in seconds. After a program has been selected further instructions may be provided. Then the actual treatment will begin.

## ***EMG-Setting***

Please follow these instructions to configure the EMG-scale in accordance with the EMG-mode set:



- **Manual EMG-Setting**  
Begin by tensing the appropriate muscle regions as strongly as possible. While doing so adjust the EMG-scale with *Button C & D* until roughly  $\frac{3}{4}$  of the scale is filled (see illustration above).
- **Semi-automatic EMG-Setting**  
Tense the appropriate muscle regions as strongly as possible. The EMG-scale will adjust automatically to the best possible range. This process takes about two seconds on average to finish.
- **Automatic EMG-Setting**  
EMG-scale will be adjusted dynamically during the exercise. As such the EMG-scale does not have to be set.

Press *Button A* to finish the EMG-setting.



## **Electrical current**

Prior to starting a training featuring electrical stimulation (e.g. ES) the intensity level of the electrical current needs to be set. The strength of the current may be adjusted during the exercise but an initial level of intensity has to be set before the training can begin.



The electrical current applied by the electrodes is displayed both as numerical values and on a logarithmical scaled graph.

To increase or decrease the intensity level of the current press *Button C & D*. A plus and minus icon are shown on the screen next to these buttons to indicate their function. Note that during an exercise the electrical current is 'locked' and cannot be increased right away. To unlock this setting press *Button D* to decreased the current first.

Check the electrodes and their contact area if the level of intensity of the electrical stimulation cannot be increased. The device is able to detect a lack of contact between electrodes and skin and will show a special warning icon on screen if this occurs.

The electrical stimulation should never be painful or uncomfortable. How the stimulation is perceived depends heavily on the current frequency and individual sensitivity to electricity.

Press *Button A* to finish the Electrical current setting.

## ***Pressure Setting and Zero Adjustment***

Depending on the sensor type a minimal initial pressure may occur when connecting a sensor with the device.

This initial pressure level is detected after the program selection is completed and compensated with a zero level adjustment. Please do not disturb or stress the pressure sensor at this point.

The pressure values displayed during the pressure setting should stabilize before the pressure sensor is inserted and the exercise continued.

Press *Button A* to finish the Pressure Setting.

## ***Exercise process***

After all necessary preparations are complete the actual exercise will begin.


The details of the training may differ according to the recommendations made by the medical specialist responsible for the treatment.

Several different animations for visualizing muscle activity are available during exercise phases which measure said muscle activity. In order to switch to the next animation simple press *Button A* during relax or tense phases of the training.


The option to switch between different animations may be disabled. In this case the device is locked with a single animation.

## ***EMG-biofeedback training***

- Follow the instructions described under *Preparations*, *Program selection* and *EMG-Setting*.
- The training consists of two alternating phases:
  - Relax




- Tense




- Follow the instructions and passively relax or actively tense the relevant muscles.

## ***Electrical stimulation training***

- Follow the instructions described under *Preparations*, *Program Selection*, *EMG-Setting* and *Electrical current*.
- The training consists of two alternating phases:
  - Break/Relax



- Electrical stimulation




- Follow the instructions and passively relax or actively tense the relevant muscles.
- The electrical stimulation phase begins with an audio announcement stating „Electricity“. Following this an electrical current will be applied to tense the musculature.
- Try to concentrate on the relevant muscles for the entire duration of the training and try to keep them as relaxed as possible.
- The intensity of the current can be changed at any point during a stimulation phase. See *Electrical current* for details.

## **Combined training (Biofeedback + electrical stimulation)**

- Follow the instructions described under *Preparations, Program selection, EMG-Setting* and *Electrical current*.
- The training consists of a recurrent series of four phases:
  - Break/Relax




- Tense



- Relax




- Electrical stimulation




- Follow the instructions and passively relax or actively tense the relevant muscles.
- The electrical stimulation phase begins with an audio announcement stating „Electricity“. Following this an electrical current will be applied to tense the musculature.
- The intensity of the current can be changed at any point during a stimulation phase. See *Electrical current* for details.

## ***EMG-triggered electrical stimulation***


- Follow the instructions described under *Preparations, Program selection, EMG-Setting and Electrical current.*
- The training consists of a recurrent series of three phases:
  - Break/Relax




- Tense



- Electrical stimulation



- Follow the instructions and passively relax or actively tense the relevant muscles
- Try to cross the red line on the EMG-scale during *Tense phases*. Once done the *Electrical stimulation phase* begins with an audio announcement stating „Electricity“. Following this an electrical current will be applied to tense the musculature.
- Note how the muscles tense up when they are being stimulated with an electrical current and attempt to replicate this during *Tense phases*.
- In case you are unable to cross the red line it will automatically lower itself after a short time if the option *Automatic* is activated. Furthermore a *time limit* may be set so that the device will switch to the *Electrical stimulation phase* on its own if the *Tense phase* once the timer runs out.
- You can change the strength of the electrical stimulation at any point (see *Electrical current*).
- If the screen displays the text *wait*, either increase the EMG-scale by pressing *Button C* or increase the trigger level threshold (red line) by pressing first *Button B* (red line blinking) and then *Button C*.





## ***Pressure measurement***

This program aims to measure the pressure experienced by the pressure sensor. The pressure decreases, when the surrounding muscles relax. All measurements are passive, meaning that no instructions for *Tense* or *Relax* phases are given during this program.

## ***Pressure measurement + 1-channel EMG***

This program combines pressure measurement with EMG-measurement on *Socket CH1*. This mode of operation is suitable for many different applications, such as simultaneous measurement of rectal pressure and tension of abdomen muscles.

To accommodate these different measurements the display was separated into two parts. The upper part of the screen displays the pressure as a line graph. The lower part shows EMG-voltages as real values on a color scale.

## ***Pressure measurement + 2-channel EMG***

This program combines pressure measurement with EMG-measurement on both *Socket CH1* and *CH2*.

The display of the program is similar to pressure measurement with only one channel EMG, except that the EMG-display is divided into two halves. The left half is shows the input of *CH1*, the right one is assigned to *CH2*.

## ***Pressure training***

This program is similar to EMG-biofeedback program. It consists of only two phases: *Relax* and *Tense*.

The *Tense* phase may be absent if its duration has been set to 0 seconds.

The EMG-scale automatically adjusts its range to fit the measurement. There are three visualizations available: line graph, bar graph and a flower animation.

## ***Finishing a therapy session***

- The training can be stopped at any time by pressing the *On/Off Button*.
- If an exercise duration has been set the device will end the training session on its own once the training duration has been reached.
- At the end of every pressure measurement a summary statistic screen will be displayed. Pressing the *On/Off Button* will turn the device off. Pressing the *Button A* will restart the device.
- Once the device has been turned off, pull the electrode cables from their sockets and remove the electrodes and sensors. Please store electrodes and sensors as described under *Maintenance*.  
Remove the pressure sensor if used from *Socket CH3*.

# 7 Maintenance

## ***Maintaining sensors***

Clean sensors after every use with running, lukewarm water. In case of tough stains, soap and washing lotion may be used in addition to water. Do not use aggressive or abrasive cleaning agents of any kind.

After cleaning the sensor dry it with a soft cloth and store it in an open plastic bag. A closed bag or storage container may amplify germ formation.


## ***Maintaining electrodes***

In case the self-adhesive electrodes no longer stick to skin, moisten the sticky side with a few water drops.

Note that before placing electrodes on skin, the relevant skin areas should be cleaned with soapy water or an alcohol wipe.

## ***Maintaining the device***

Use a dry, soft cloth to clean the device of light stains. Do not use aggressive or abrasive cleaning agents of any kind. In case of tougher stains use soapy water or washing lotion in addition to a soft, moist cloth. Take care that no water enters any of the sockets or battery compartment. Dry the device afterwards, again with a soft cloth.

	<b>CAUTION</b>
<b>Hazard due to defective device</b>	
→ Never clean the device under running water.	
→ Keep sockets and cable jacks dry at all times.	

## Replacing batteries

Under ideal circumstances one set of batteries will last a Conpreva<sub>2</sub> for approximately 10 hours of use.

When the battery charge is getting too low the device will announce this with a double beep warning before turning off.

By default the device is delivered with one set of batteries. To replace the batteries unlock and open the *Battery compartment* on the back of the device.

**Pay close attention to the polarity of the new batteries as you place them within the device.** See picture below for correction battery orientation.



**Once all batteries have been replaced: close the battery compartment again.**

Rechargeable batteries are available as optional accessories. These can be recharged with an external battery charger. Removing and placing them back in the device is done in the same way as with regular batteries. Please only recharge these batteries once they are completely empty, as indicated by the device turning off with two beep noises.

Recharging batteries when they still carry much charge can damage them and make it necessary to replace them early on.

## Safety checks

The Conpreva<sub>2</sub> has been manufactured in accordance with the Medical Devices Act and accordingly does not require any safety checks. For details of the applicable standard in your country please contact the local distributor or the manufacturer of this medical device.

## Troubleshooting

See the table below if the device shows any errant behavior. Please note that only specially trained persons may conduct repairs on the device and its accessories.

Error	Cause	Solution
No measurable EMG-signal	No connection between sensor cable and <i>sensor socket</i>	Check that the cable plug is correctly connected with the <i>sensor socket CH1</i>
	Insufficient contact between skin and electrode	Use new electrodes or use contact gel on the sticky side of the electrodes
	EMG-Range is set too high	Select a smaller EMG-Range (see <i>EMG-Setting</i> )
EMG-Signal overflow (scale	No EMG-signal – see above	See above

completely filled)	EMG-Range is too small	Select a larger EMG-Range (see <i>EMG-Setting</i> )
	Strong electromagnetic fields are causing interference	Change location to somewhere clear of such interference
Device lasts <2h with fully charged batteries	Batteries are defective	Replace batteries
Device does not turn on	Batteries are fully discharged	Replace/charge batteries
During the stimulation phase of a training session a yellow warning sign appears	Sensor cable is not plugged in correctly.	Check that the cable plug is connected with the <i>sensor socket CH1</i>
	Insufficient contact between skin and electrode	Use new electrodes or use contact gel on the sticky side of the electrodes
	Device has reached its limits (very high skin resistance)	Use contact gel when attaching electrodes
Intensity of the electrical current does not increase	The device locks the electrical current during training to prevent its accidental increase.	Press <i>Button D</i> to unlock the electrical current and <i>Button C</i> to increase it manually
No electrical current is being applied during the training	Wrong mode of operation for the training	Mode of operation does not include electrical stimulation
	Current intensity is too low	Increase the electrical current

## Storage and transport



### CAUTION

#### Danger to user due defective device



→ Medical device may only be stored at an ambient temperature between 5°C and 35°C.



→ Do not expose the device to water.



→ Protect the device from strong sunlight or heat sources.

→ Do not use in dusty environments.



→ Medical device may only be stored at a relative humidity between 30-75%.



→ Medical device may only be stored at a relative pressure between 700-1060 hPa.

## **8 Environmental protection & disposal**

This medical device conforms to the requirements of guideline 2011/65/EU of the European Parliament and Council (dated June, 8, 2011 and amendments) on the restriction of certain hazardous substances in electrical and electronic equipment.

This device may not be disassembled or by disposed as household waste once the treatment period is over. Please return the device to its manufacturer or supplier for proper disposal.

Under proper care the Conpreva<sub>2</sub> has a product life cycle of up to 10 years. This does not exclude that repairs may be necessary during that timeframe.

### **Batteries and accumulators**

Please dispose of batteries according to your local or national legal provisions (within the EU: guideline 2006/66/EG (dated September 6, 2006) of the European Parliament and Council).



## 9 Technical information

### *Name of the medical device*

HPS12

### *Model*

Conpreva<sub>2</sub>

### *Display*

LCD color screen

Resolution: 320x240 pixel

### *Connection ports*

- 1x 3-pole socket on the front side to connect electrode cables or sensors for EMG-biofeedback or electrical stimulation (CH1)
- 1x 3-pole socket on the front side to connect electrode cables or sensors for EMG-biofeedback (CH2)
- 1x 4-pole socket for PC connection cable (USB, 5000V checked/ optical decoupled)
- 1 x Luer-Lock adapter, female (CH3) for connecting a pressure sensor at the front side of the device
- 1 x 2-pole socket on left side of device for connecting an optional power supply unit (6V  $\overline{\text{---}}$  / 1A; 1.3mm Stift;  $\varnothing$  3.7mm max). The power supply unit must be conform with **EN60601-1:2006** and **EN60601-1-2:2007**.

## ***Modes of operation***

- EMG-biofeedback
- Electrical stimulation
- Combined training (Biofeedback + electrical stimulation)
- EMG-triggered electrical stimulation
- Pressure measurement
- Pressure training
- Pressure measurement + 1-channel EMG
- Pressure measurement + 2-channel EMG


## ***EMG-Range***

- 0 – 300  $\mu\text{V}$  with a resolution of 0.01 $\mu\text{V}$  ... 0.25  $\mu\text{V}$
- Displaying of measurements both as numerical values and animations
- Sound rendering of EMG-measurements
- EMG-Range can be set manually, semi-automatic or fully automatic


## ***Output***

Stimulation current	mono/biphasic
Maximum current	250mA at 500 Ohm at 27Hz
Root mean square current	9.9mA at 1000 Ohm at 9.9V
Pulse width	Adjustable: 40 - 500 $\mu\text{s}$
Frequency	Adjustable: - 2 Hz - 200 Hz LF - 25 kHz MF
Max. peak at 500 Ohm	$\pm 160\text{V}$ at 5Hz
Max. peak at open output	$\pm 170\text{ V}$
Max. energy of a single pulse with a load impedance of 500 Ohm	2 mJ

## ***Pulse pattern***



## ***Pulse group pattern (medium frequency)***



## ***Exercise phases***

- Stimulation phase: 2...255s
- Tense phase: 0...255s
- Break phase: 0...255s
- Relax phase: 1...255s
- Exercise duration: 1...60min & continuous

## ***Power supply***

- 4x mignon AA batteries  
*or*
- 4x NiMH rechargeable batteries (2500mAh, 1.2V)  
*or*
- Power supply unit (6V  $\text{---}$  / 1A; 1.3mm plug;  $\varnothing$  3.7mm max). The power supply unit must conform with **EN60601-1:2006** and **EN60601-1-2:2007**.

## ***Safety features***

- Electrical current cannot be increases outside of stimulation phases
- Gentle rise and fall of the electrical stimulation
- Warning upon electrodes having too little contact with skin (electrical current is limited and adjusted)
- Electrical current is locked by default and needs to be unlocked before it can be increased (lock reengages after 8 seconds without change)
- Protection class II, IP 20

## **Attention**

Please mind the storage conditions set for accessories.

## Operating conditions



### CAUTION



#### **Danger to user due to defective device**

→ Medical device may only be used at an ambient temperature between 5°C and 35°C.



→ Do not expose the device to water during therapy sessions.



→ Protect the device from strong sunlight or heat sources during therapy sessions.

→ Do not use in dusty environments.



→ Medical device may only be operated at a relative humidity between 30-75%.



→ Medical device may only be operated at a relative pressure between 700-1060 hPa.

---

## **Dimensions and weight**

- L x W x H: 95x150x30 mm
- Weight (excl. batteries): 250 g

## **Storage of exercise data**

Exercises are recorded and their data stored on the device. This includes the date, time and mode of operation of a training exercise as well as each individual EMG-value measured. Thus it is possible to check if the training schedule has been kept or not. Each patient is assigned only training sessions as conducted by them. This permits an objective analysis of the treatment process by the medical specialist responsible. The device is capable of managing up to six different medical specialist, whom may use the device, each of which can be assigned up to six patients. Every patient may have up to six different training programs. Reading the stored data from the device should only be done by authorized persons (see manual for HPS-Software).

## ***Application recommendation***

EMG-biofeedback treatment with the Conpreva<sub>2</sub> allows for selective training of muscles to enhance their coordination, endurance and strength.

Optimal training results can be reached with strength training by using semi-automatic EMG-Setting. Simply tense the relevant muscles as hard as possible for about 5 seconds while the EMG-Range is being set.

# **ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES**

Cable length of **K1 (MED) electrode cable**

	Cable length	Recorded length
K1med electrode cable	1.5 m	1.5 m

The use of longer cables can lead to increased electromagnetic emissions or reduced immunity to interference. Only use original manufacturer accessories.

## **Guidance and manufacturers declaration - electromagnetic emissions (Tab. 201 in accordance to DIN EN 60601-1-2)**

The device is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment guidance
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The device is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2		Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3		Not applicable

**Guidance and manufacturers declaration — electromagnetic resistance  
(Tab. 202 in accordance to DIN EN 60601-1-2)**


The device is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601- Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Magnetic field with a supply frequency (50/60 Hz) IEC 61000-4-8	3 A/m	Not applicable because magnetic sensitive	



**Guidance and manufacturers declaration — electromagnetic resistance  
(Tab. 204 in accordance to DIN EN 60601-1-2)**

The device is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.

Immunity level	IEC 60601 - Test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile radio devices should not be used any closer to the device or any part of it, including cables, than the recommended minimum distance calculated from the equation applicable to the frequency of the transmitter. Recommended distance:
Conducted HF-disturbance IEC 61000-4-6	3 Veff 150 KHz - 80 MHz	$V_1 = 10 \cdot V_{emk} @$ 150 KHz to 80 MHz	$d = (3.5/V_1) \sqrt{P}$
Radiated HF-disturbance IEC 61000-4-3	3 V/m 80 MHz - 2.5 GHz	$E_1 = 10 \text{ V/m} @$ 80 MHz to 2.5 GHz	$d = (3.5/E_1) \sqrt{P}$ for 80MHz to 800 MHz
			$d = (7/E_1) \sqrt{P}$ for 800 MHz to 2.5 GHz  Where P is the maximum output rating of the transmitter in watt (W) according to the transmitter manufacturer and d is the recommended separation distance in meter (m). Field strengths from fixed radio transmitters, as determined by an electromagnetic survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80MHz and 800 MHz, the higher frequency range applies.  
 NOTE 2: These guidelines may not apply in all situations.  
 The spread of electromagnetic fields is affected by absorption and reflection from structures, objects and people.

a) The field strength of fixed transmitters, such as cellular/cordless telephone base stations, mobile land radios, amateur radios, AM and FM radio broadcasts and TV broadcasts cannot be predicted with sufficient accuracy. Investigations are required to assess the electromagnetic environment at a location with any precision. If the measured field strength at the place the device is used exceeds the compliance levels described above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b) Beyond the frequency range of 150kHz to 80MHz, field strength should be less than 10 V/M.

**Recommended distances between portable and mobile high-frequency communication equipment and the device (Tab. 206 in accordance to DIN EN 60601-1-2)**

The device is intended for use in an electromagnetic environment with limited high frequency disturbances. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radio communication equipment (i.e. transmitters) and the device as recommended below, in accordance with the maximum output of the communications equipment.

	Protective distance depending on the frequency of the transmitter in meter		
Rated maximum output of transmitter in W (Watt)	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d=(3.5/V_1)\sqrt{P}$	$d=(3.5/E_1)\sqrt{P}$	$d=(7/E_1)\sqrt{P}$
0.01	0.12	0.04	0.08
0.1	0.35	0.11	0.22
1	1.11	0.35	0.75
10	3.51	1.11	2.22
100	11.90	12.00	22.00
For transmitters rated at a maximum output not listed above, the recommended protective distance in meters (m) can be estimated by applying the equation to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. The spread of electromagnetic fields is affected by absorption and reflection from structures, objects and people.			

## 10 Warranty

The original manufacturer provides statutory warranty for the device and all accessories produced by this same manufacturer from the date of the original consumer purchase.

If the product fails to function properly during its warranty period due to a defective part or faulty workmanship, the manufacturer or selling dealer will replace or repair the product.

All repairs must be performed by a service center authorized by the manufacturer. Any attempts to modify or repair the device performed by unauthorized persons will void this warranty.

Any malfunction or failure of the device caused by misuse of the product, including but not limited to, failure to provide reasonable and necessary maintenance, use inconsistent with the product's user manual or deviating from the instructions given by medical specialists is also not covered by warranty.

If problems with the device or its accessories arise, please contact the manufacturer:

Haynl-Elektronik GmbH  
Magdeburger Str. 117a  
39218 Schönebeck  
Germany  
Tel.: +49 (3928) / 69414  
Fax: +49 (3928) / 76222  
[www.haynl.com](http://www.haynl.com)  
[info@haynl.com](mailto:info@haynl.com)


# 11 Advanced information for physicians and therapist



## CAUTION

### Danger to users due to incorrect therapy parameters

→ Only authorized persons may access advanced information and therapy settings

- The Conpreva<sub>2</sub> offers extensive options by which an optimal and personalized therapy plan can be designed for each individual patient.
- At the core of these options sits the main menu (see image below).
- To access the main menu press *Button D* while the device is booting. A gear icon is displayed next to *Button D* on the screen during the boot sequence to indicate this.
- Use the touchscreen to navigate the main menu, access sub-menus and configure options.
- Press *Button D* to leave a submenu or return to a previous screen. On screen an arrow icon  is shown next to this button (see image). The icon itself has no touch function.
- To exit the main menu press *Button D* while the first screen of the main menu is shown (see image below).
- It is necessary to enter an access code if the option code request is enabled. Enter the code via the touchscreen when a number block appears on screen.
- The access code is enclosed as a leaflet at the end of this manual. Please keep the access code secret so that patients do not unwittingly change therapy parameters.
- An invalid code will cause the device to restart.

The main menu is divided into three major sub-menus:

### **Therapy Settings**

This sub-menu contains all parameters needed to fully configure therapy programs, such as mode of operation, duration of each individual trainings phase and properties of the electrical stimulation.

Treatment patters for common disorders or syndromes are predefined. Based on their diagnosis and up to date medical knowledge optimal treatment parameters can be devised for each individual patient.

### **System Settings**

The options in this sub-menu concern more general aspects of the device. This includes the sound volume and animations used during training sessions as well as a number of parameters for operating the device.

### **Administration**


This menu all therapists using the device, their assigned patients and stored exercise data of these patients can be managed.



## Therapy Settings

This sub-menu provides extensive functions to configure training programs so that they fit the specific treatment of each individual patient who will use the device.

Initially a list of the programs assigned to the currently selected patient will be provided. Which patient is selected and the therapist they belong to is displayed in the upper right corner of the screen. *Buttons B & C* will switch forward through the list of patients and therapists. Once a program is selected a list parameters to be configured will be displayed (see image below). These parameters can be configured:



### Mode of operation

This is the central property of any training program and should be set first. It defines whether a program will feature EMG-biofeedback (EMG), electrical stimulation (ES), a combination of these two (EMG+ES) or EMG-triggered electrical stimulation (Triggered).

### Current settings

This submenu provides access to all parameters of the electrical stimulation. This includes *current frequency*, *pulse width* and *waveform*.

## Phases

Each exercise consists of multiple repeating phases which includes phases for relaxing, tensing muscles, electrical stimulation and a follow break period. The duration of these training phases can be individually set in this sub-menu.

The duration of the entire training session with this program can also be set here. The maximum duration is one hour, beyond that training can be set to be continuous so that it only ends when the device is turned off.

## Program name

This property defines the name assigned to the training program. It will be displayed in the *program selection* when patients have to choose which program they want to start. The name is also used in the list of programs shown in the *therapy settings* sub-menu.

The name is entered via the on-screen keyboard displayed after accessing this property.

The keys labeled [  $\wedge$  ] and [123] on the left side of the keyboard are used to switch to capital letters and number block.


To confirm the input simply return to the previous menu.

## Electrode positions

This property defines if and which images for placing sensors or electrodes are shown to patients before the start of an exercise.

## System Settings

This sub-menu includes the numerous settings defining the general behavior of the device:



### Multimedia

In this menu the sounds played during training sessions, namely *EMG-sound* and *voice output*, can be enabled. The *volume* of these sounds is set via a separated slider.

### Therapy Parameters

Some therapy settings apply across multiple training programs and as such are set globally.

### EMG-Setting

The *EMG-Setting* defines the method by which the EMG-Range is calibrated before the start of a training session. Available options are *Manual*, *Semi-automatic* and *Automatic*. A definition of these different modes can be found under section *Training* → *EMG-Setting*.

### EMG Maximum

This option sets an upper limit for the EMG-Range.

### Trigger Settings

The *Trigger Settings* define the behavior of the trigger threshold (red line visible during tense phase) in programs featuring EMG-triggered stimulation.



Please note that when both *Automatic Trigger* and *Time Limited Trigger* are disabled, the Tense phase may continue indefinitely unless the trigger threshold is reached. This may cause a deadlock.

***Automatic Trigger***

When this option is enabled the electrical stimulation will commence even if the patient did not reach the trigger threshold.

***Time Limited Trigger***

With this option active the trigger will gradually lower itself each time a pre-defined time period has passed without the trigger threshold being crossed.

***Time Limit***

This time limit defines the maximum time available to a patient for crossing the trigger threshold during tense phases (up to 200s). *Automatic Trigger* has to be enabled for this to have an effect.

***Update rate***

This value defines the time period (max. 50s) which has to pass before the trigger threshold is updated to a lower value. This option has no effect without *Time Limited Trigger* enabled.

## **Device Parameters**

This menu contains a number of options which define the general behavior of the device.

### **Patient Selection**

If this option is enabled, the device will show a patient selection screen prior to the program selection screen. On that screen the user can choose the patient to whom the training data is to be assigned. Note that patients are grouped by their therapist.

In case that this option is disabled all training data is assigned to the first patient of the first therapist instead. The device will also skip the patient selection screen and begin with the program selection screen instead (if more than one training program is active for the patient).

### **Access Code**

This option enables or disables the access code request prior to accessing the main menu.

Please note that unauthorized persons may change therapy and device settings if no access code is required to enter the main menu.

### **Button Feedback**

This option enables or disables the feedback sound when pressing a button or the touchscreen.

### **Lock Animation**

If this option is enabled, the patient will not be able to change the animation shown during training.

### **EMG-Display**

This option enables or disables the numerical display of EMG-value during training.

### ***Numeric Display***

This option enables a numerical display for EMG-values during exercises.

### **Show EMG after current**

This option enables the EMG-scale display during the Break phase after a stimulation phase.

### **Pressure Unit**

This menu sets the default measurement unit used for displaying pressure values. Available units are *mbar*, *mmHg* and *psi*.

### **Maintenance**

This menu features a number of technical options and functions for managing the data stored on the device.

### **Extras**

This submenu contains settings such as the *Baud Rate* (data transmission speed) for PC-connections, *Screen Brightness* and toggling *Bluetooth* connectivity.

### **Delete Exercise Data**

With this function all exercise data stored on the device can be deleted at once. This requires a manual confirmation to prevent unintended loss of data.

Please note that with this action all exercise data will be lost with no possibility of recovery.

### **Factory Reset**

This function will reset the device completely to factory defaults. As such all exercise data, training programs, patient profiles therapist profiles and settings are lost. This includes all therapy and device settings.

Please note that this process is non-reversible and all data, programs and settings will be is lost with no possibility of recovery.

### **Language**

Here users may change the language used both during training sessions as well as within the menu.

## Administration

In this menu both therapists and patient profiles are managed. It also provides detailed statistic data for every exercise stored. Initially the screen will show a list of therapists (see image below).



## Therapists

Press any label with a therapist's name to access their profile. Here their *name*, *surname* and *department* can be changed. To do so, press the button with a pencil icon next to the property.

## Patients

To view the patients assigned to a therapist, press the button with the group icon next to the relevant therapist name. The screen will show a list of patients who are assigned to the therapist chosen. Accessing and changing a patients details is done the same way as with therapists.

Additionally, if there is any exercise data assigned to a patient, a button will be displayed next to their name in the list of patients. Pressing this button will display a list of this patient's recorded exercises. Pressing on the name of any such record will display details and statistical data of the exercise in turn.

Any exercise can be deleted by pressing the button with the waste basket icon. Scroll through the list of recorded exercises by dragging your finger over the vertical scroll bar.







**Manufacturer:**

*Haynl-Elektronik GmbH  
Magdeburger Straße 117a  
D-39218 Schönebeck  
Germany*

