CE 2797

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GB

Table of contents

Pictograms **40**
General information **41**
Important note **42**
... regarding these operating instructions 42
... regarding the intended purpose/operating instructions 42
... regarding liability and warranty 42
Clinical usage 43
User qualification 43
Safety instructions 43
Residual risks 47
Contraindications/Application restrictions 47
Device description **49**
TENS pain therapy with SaneoTENS 49
Indications 49
Program overview 50
Program overview by indication/body region 52
Scope of delivery 54
Accessories 54
Device images 56
Possible displays 58
Operating SaneoTENS **60**
Prepare application 60
Place the electrodes 60
Switching SaneoTENS on and off 61
Set the program 61
Select stimulation intensity 62
Programming freely available programs 63
End treatment, cleaning 65
Battery change 66
Technical data and support **67**
Technical specifications 67
Labelling 68
Disposal / Maintenance 68
Fault and problem handling 69
Notes **71**

Pictograms

GB

	Make sure that you read through the operating instructions before first use and observe all instructions.
	Earth-free applied part **B**ody **F**loating type. The product affords protection against electric shock at permissible leakage current.
	This product bears the recycling symbol according to the EU Directive 2002/96/EC. It is not permitted to dispose of the product in the household/ residual waste. The contact partner responsible for the disposal is the dealer/manufacturer.
CE 0000	With the CE marking, manufacturers indicate that their product complies with the requirements of a European directive or an EU regulation. With the CE marking, the medical device manufacturer expresses the conformity of the product with EU Regulation 93/42/EEC (MDD) or Regulation (EU) 2017/745 (MDR). If a four-digit identification number is placed after the CE marking, this indicates the involvement of a notified body.
	Manufacturer
	Year of manufacture
SN	Serial number of the device
LOT	Batch number
	Usable until
MD	Medical device

CAUTION	Identifies a possibly harmful situation. If it is not avoided, property damage may be caused.
⚠ **CAUTION**	Indicates a potentially dangerous situation. If not avoided, injuries may result.
	User information: gives important advice, e. g. regarding the configuration.

General information

Dear customer,

many thanks for your trust and for deciding on this high-quality product from

tic Medizintechnik GmbH & Co. KG.

These operating instructions will guide you through the handling and application options of the SaneoTENS therapy system.

Please always keep these operating instructions with the device. Instructions for use suitable for blind and visually impaired people (e.g. in MS Word format [read aloud function]) can be provided upon request.

This device does not generally replace a visit to the doctor!

As the use of the therapy device has shown good results, despite significant pain relief, this should not be equated with a cure for the cause of the pain. We recommend that you consult your doctor regularly.

Device properties

- 19 therapy programs (3 of which are user-definable)
- Timer function
- Large LC display with complete and clear display of functions and parameters
- Automatic switch off after the end of treatment
- Easy use
- Battery monitoring
- Automatic electrode detection
- Lock function of the current strength
- Belt clip for more mobility

Important note

GB

... regarding these operating instructions

Please read through these operating instructions before first use.

These operating instructions are part of the device and must be available to the user at all times. Extracts are not permitted.

Please be sure to observe the instructions for intended use (page 42), the safety instructions (page 43 et seq.) and the contraindications (page 47).

... regarding the intended purpose/operating instructions

The SaneoTENS is a home therapy device used for electrical nerve stimulation through the skin. Please use it only for the treatment of the indications described in the **Indications** chapter and in the program overview (page 49 et seq.).

Any use that is not in the operating instructions or described as not permitted is unintended use.

Only accessories authorised by the manufacturer may be used (e.g. TENS body surface electrodes).

... regarding liability and warranty

From the date of delivery, the SaneoTENS therapy device and accessories (excluding consumables) are covered by a 24-month warranty. During this period the device/accessories will be repaired or replaced free of charge by the manufacturer.

Unintended use, failure to follow these instructions and unauthorised modifications to the system exclude the manufacturer's liability for resulting property damage and personal injury. This voids the warranty.

Important note

The manufacturer reserves the right to make technical changes to improve the therapy device, its functionality and accessories.

Clinical usage

The clinical benefit of SaneoTENS lies in the nerve stimulation triggered by electrical impulses, which allows the user to achieve pain relief in the affected area of the body, depending on the program selected and tailored to their specific indication.

User qualification

Before using the SaneoTENS for the first time, it is recommended that the user be instructed by medically qualified personnel.

It is assumed that

- the patient is at least 18 years of age
- andhas the cognitive skills to administer the treatment him/herself and to read and understand the instructions for use.

Safety instructions

The SaneoTENS is a class IIa medical device.

Only use the SaneoTENS

- with original accessories.
- protect against water or other liquids and do not use at a relative humidity of over 90 %.
- do not use, if it is faulty or damaged (e.g. from a fall).
- keep it in the packaging and protect it against damage and dirt.
- do not use it while driving or operating machinery (e.g. driving car) or while sleeping.

Important note

- do not operate it with the battery compartment open.
- do not use it in water.

GB

Please observe the operating and storage conditions as described in the chapter **Technical data**.

CAUTION **European standard**

This device complies with European standards EN 60601-1 and EN 60601-1-2 (compliance with IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-6) and is subject to special precautions regarding electromagnetic compatibility.

CAUTION **Electromagnetic compatibility**

Using the device in the presence of electromagnetic interference can lead to error messages or a failure of the display/device.

CAUTION **Do not stack the device!**

Use of this device in close proximity to other devices or in a stacked configuration should be avoided as this could result in improper operation. If use in the manner described above is still necessary, this device and the other devices should be monitored to ensure that they are functioning properly.

CAUTION **Electric power source**

Remove the batteries from the device if you will not be using it for a long period of time. Leaking batteries can damage the product.

Important note

CAUTION **System malfunction due to broken electrode cable!**

A damaged electrode cable can lead to system malfunctions (e.g. due to improper use and extreme bending or sudden pulling).

➔ Do not bend the electrode cables.

➔ Do not pull on the cables.

CAUTION **Interactions with other devices/systems!**

The output values can be influenced by interactions with other (medical) devices. Do not use the SaneoTENS near:

- Microwaves or short-wave devices.
- Portable HF communication devices (e.g. mobile phones)

➔ During treatment please keep your mobile phone at least three metres away.

CAUTION **Only use certified accessories!**

The use of accessories other than those specified or provided by the manufacturer of this device may result in increased electromagnetic interference or reduced electromagnetic immunity of the device and may lead to improper operation.

CAUTION **Application limitation**

Do not apply the electrical stimulation treatment over or through the head, directly on the eyes, covering the mouth, on the front of the neck or on the chest and upper back, or crossing the heart.

Important note

GB

⚠ CAUTION **Uncontrolled current output!**

Touching metallic parts in the battery compartment can lead to uncontrolled discharges during application. Only operate the device with the battery compartment closed!

⚠ CAUTION **Use under medical observation**

The patient/user may only perform the treatment under medical observation, if:

- mental limitations or cognitive differences are known in the patient/user:
- feverish diseases, local inflammations, thrombosis, inflammatory skin changes are present in the area of application.

⚠ CAUTION **Ventricular fibrillation!**

Application to the chest and near the carotid artery may increase the risk of ventricular fibrillation!

⚠ CAUTION **Danger of strangulation from electrode cables!**

Keep the device and the electrode cables in an inaccessible place for children and unauthorised persons and do not use it in your sleep.

⚠ CAUTION **Electrochemical burns or irritation!**

The SaneoTENS is not to be used with simultaneous connection to HF surgery or electrocoagulation devices. Non-observance could lead to electrochemical burns / irritation.

Important note

Residual risks

The residual side effects identified in the literature on electrostimulation devices after implementation of all risk-reducing measures are described below.

Product: The following remaining side effects have been reported in the clinical literature for equivalent products and may be associated with SaneoTENS itself:

Difficulties in handling the device

Method: The following residual side effects have been reported in the clinical literature for treatment with electrostimulation devices and may be associated with the method itself:

Redness/swelling at the site of surface electrode placement, pain, discomfort, pressure sensation, hypersensitivity or irritation

Contraindications/Application restrictions

Medical conditions, under which the SaneoTENS device may **not** be used:

- During pregnancy never use the device in the abdomen region. Medically effective current can lead to a miscarriage of an embryo!
- While operating machinery (e.g. driving a car) or sleeping.
- With metallic implants in the treatment area.
- In the case of skin and mucosa defects (injuries) and acute inflammation of the mucosa, skin or subcutaneous tissue in the treatment area.
- In the case of cardiac arrhythmias or heart muscle diseases.

Important note

GB

- If mental limitations or cognitive deficits are known.
- Never use in the heart area or near the carotid sinus (carotid artery).
- The device is not suitable for persons with pacemakers and/or active implants (defibrillators)!

- Medically effective current promotes blood circulation. This can lead to reddening of the skin in the treatment area. Treatment must be interrupted and/or carried out under medical supervision if:
 - the current sensitivity changes drastically, even though the stimulation intensity was not changed.
 - unknown inflammation or bleeding in the treatment area is noticed before the start of treatment.
 - increasingly larger skin changes occur during treatment.
- Allergic reactions to tangible materials possible (for materials, please refer to the list of scope of delivery in the device description of these instructions)

Device description

TENS pain therapy with SaneoTENS

TENS (**T**ranscutaneous **E**lectrical **N**erve **S**timulation) has been recognised and proven in medicine for over 50 years. The SaneoTENS is a TENS therapy system for the effective treatment of acute and chronic pain.

During TENS therapy, a medically effective and harmless current is conducted through the skin into the tissue via two or four electrodes.

Two frequency ranges are distinguished in the treatment programs.

- **TENS therapy (2 Hz to 20 Hz)**

 ... is ideally suited for the treatment of chronic pain because it can stimulate the release of pain-relieving endogenous substances (e.g. endorphins, dynorphins). The body's own painkillers are also effective beyond the duration of the treatment.

- **TENS therapy (50 Hz to 100 Hz)**

 ... is particularly suitable for the treatment of acute pain. It can block or reduce pain transmission to the brain.

Indications

Pain conditions that cannot be influenced by the cause, e.g.:

- Sympathetic pain syndrome
- Post-herpetic neuralgia
- Deafferentation stump phantom pain
- Central pain
- Chronic pain after CNS lesion
- Secondary central pain, chronic radicular and secondary chronic pain processes

Program overview

GB

The following table lists all **programs** with their respective **areas of application**, **mode of action**, **settings** & **duration**.

Programs	Areas of application	Mode of action	Frequency (Hz)	Impulse width (µs)	Increase (sec)	Plateau (sec)	Decrease (sec)	Pause (sec)	Duration (min)
1	Acute pain, nerve-related (neuropathic) or chronic pain	Continuous stimulation, channels A and B identical	100	200					30
2	Acute pain, nerve-related (neuropathic) or chronic pain	Cyclic stimulation, channels A and B identical, 3 sec 100 Hz, 3 sec 3 Hz	100/3	180					30
3	Alternative to programs 1 and 2, especially for chronic pain	Continuous stimulation, Channel A Channel B **(apply channels A and B crosswise)**	 100 4	 200 200					30
4	Alternative to programs 1, 2 and 3 for the prevention of therapy resistance	Burst stimulation, channels A and B identical, 7 pulses per burst group, pulse frequency 2 Hz, burst frequency 100 Hz	100	250					30
5	Alternative to programs 1, 2 and 3 for the prevention of therapy resistance	Frequency modulation, channel A and B identical	8-100	250	2	2	2	2	30
6	Muscle pain	Continuous stimulation, channels A and B identical, 10 minutes at 100 Hz, then 10 minutes at 4 Hz, then 10 minutes at 100 Hz.	100/4	150@ 100Hz 200@ 4Hz					30
7	Tensioning of muscles	Cyclic stimulation, Channel A and B identical, 30 sec 4 Hz, 30 sec 6 Hz, 30 sec 8 Hz, 20 sec 4 Hz, 20 sec 6 Hz, 20 sec 8 Hz, 2 second pause between each phase	4/6/8	250					30

Programs	Areas of application	Mode of action	Frequency (Hz)	Impulse width (µs)	Increase (sec)	Plateau (sec)	Decrease (sec)	Pause (sec)	Duration (min)
8	Tensioning of muscles	Amplitude modulation, channel A and B alternating	75	150					30
9	Tensioning of muscles	Cyclic stimulation, channels A and B identical, 5 sec 40 Hz, 10 sec 4 Hz, 5 sec 50 Hz	40/4/50	250					30
10	Spinal syndromes (cervical spine, thoracic spine, lumbar spine)	Continuous stimulation, channels A and B identical	110	180					30
11	Spinal syndromes (cervical spine, thoracic spine, lumbar spine)	Continuous stimulation, channels A and B identical	80	150					30
12	Knee, joint or tendon problems, phantom pain	Continuous stimulation, channels A and B identical	40	200					30
13	Rheumatic joint problems	Frequency modulation, channel A and B identical	2-8	250	2	2	2	2	30
14	Neuralgia, shingles	Continuous stimulation, channels A and B identical	10 2	200 200					30
15	Menstrual problems	Continuous stimulation, channels A and B identical	8	200					30
16	Circulatory disorders (Kaade system)	Continuous stimulation, channels A and B identical	2	250					30
u1	Free, see chapter **Programming freely available programs** Default	Base program: continuous TENS stimulation, configurable frequency and therapy time	2-120 20	200 200					5-99 30
u2	Free, see chapter **Programming freely available programs**	Base program Frequency modulation, frequency, increase, hold, decrease and therapy time configurable	2-100	100-250	0.5-5.0	1-10	0.5-5.0	1-10	5-99
u3	Free, see chapter **Programming freely available programs**	Base program Burst stimulation, Channels A and B identical, 7 pulses per burst group, Pulse frequency 2 Hz, Burst group interval and therapy time configurable	100	250				0.1-5.0	5-99

GB

Program overview by indication/body region

Programs 1–5 can generally be used to treat almost all types of pain. The following table also lists the most common indications or the areas of the body where pain most frequently occurs and provides recommendations on which additional programs you can use to treat them and where to place the electrodes on the body. **Please also see the fold-out cover at the end of this manual.**

To achieve maximum success from stimulation, it is important to place the electrodes correctly. The optimum positioning varies from person to person. So feel free to experiment a little until you feel the most pleasant and effective stimulation.

To ensure good conductivity, you should first wash and dry the affected area of skin before applying the electrodes without overlapping. The skin must also be healthy.

Indication / Body region	Program	Electrode placement
Achilles tendon	12	40
Abdominal muscles	6	19, 20, 21
Circulatory disorders (Kaade system)	16	28
Elbows	12	24
Gluteal muscles	6	29
Shingles	14	*
Hips	12	30, 31
Sciatic nerve	14	15, 16, 17
Knee joint	12	35

Device description

Indication / Body region	Program	Electrode placement
Menstrual problems	15	34
Neck and back	6	2, 4, 5, 6, 7, 8, 9, 11
Upper arm	6	22, 23
Thighs	6	32, 33
Rheumatic joint problems	13	*
Shoulder	6, 12	10, 12, 39
Tension headache	7, 8	2
Ankle joint	12	38
Trigeminal neuralgia	14	1a, 1b, 1c
Forearm (tennis elbow)	6	25, 26, 27
Lower leg	6	36, 37
Tension	7, 8, 9	*
Spine (cervical spine, thoracic spine, lumbar spine)	10, 11	5, 11, 14, 15, 16, 18

* Place the electrodes on or around the painful area.

Device description

GB

Scope of delivery

Art. no.	Description	Material
MG5005	1 x SaneoTENS control unit 1 x Operating instructions	ABS, silicone, polycarbonate Paper
	1 x Block battery (9 V - alkaline)	–
	1 x TENS electrodes 50 x 50 mm (4 pcs)	TPE (thermoplastic elastomer)
ME5004	2 x Cable for connection to body surface electrodes, 2-pin, 2 mm plug	Copper wire, PVC


Accessories

Art. no.	TENS body surface electrodes:	Material
	32 mm round (4 pieces)	TPE
	50 x 50 mm (4 pcs)	
	50 x 90 mm (4 pcs)	
Art. no.	**Textile electrodes**	**Material**
	Elbow electrode	Nylon, Lycra, silver thread
	Knee electrode	
	Stimulation sock (1 piece) Size 33 x 9 cm ± 1 cm	
	Stimulation glove (1 piece) Size 24 x 11,5 cm ± 1 cm	

Device description

Accessories illustration

TENS body surface electrodes:


Device description

Device images

GB

Overall view

1 Connecting socket electrode cable, channel A

2 Connecting socket electrode cable, channel B

3 ▲: Increase channel A/B intensity

4 ▼: Decrease channel A/B intensity


5 Menu/Confirmation button (select program)

6 On/off button

Device description

Display

All displays are listed in the chapter **Possible displays**, on page 59.

1	Battery charge state	**7**	Stimulation intensity channel B
2	Frequency mode is active	**8**	Stimulation intensity channel A
3	Current stimulation phase (flashing)	**9**	Burst-mode is active
4	Remaining treatment time in minutes:seconds	**10**	Invalid entry with active lock
5	Lock mechanism is active	**11**	Program set
6	The therapy session was interrupted due to insufficient contact of the electrodes with the skin		

Device description

Possible displays

Display	Description
Pxx	Treatment program (active program)
	**Locking**. After 20 seconds without entry, the keypad will be protected against accidental entry. Unlocking only by pressing the button ▼.
	**Invalid input / invalid button press Lock** remains active
	**Therapy session interrupted due to insufficient patient contact**. See chapter on **Fault and problem handling** (page 69), points 2 and 3
	**Treatment time** per session in minutes **(min)** and seconds **(sec)**
mA	**Stimulation intensity** from 0 mA to 60 mA
Burst	Mode against habituation tendencies
FM	Frequency modulation
◢■◣	The current stimulation phase or rest phase will be displayed by flashing. Image from the left: _Pause \| ◢Increase \| ■ Plateau \| ◣ Decrease \| Pause_
Units of the values shown	
Hz	Frequency (Impulses per second)
µs	Impulse duration (in microseconds)
sec	Stimulation or rest phase (in seconds)
min	Treatment duration (in minutes)

Possible displays

Display	Description
Battery charge state	
	The battery is full, device is ready for operation.
	The device is still ready for operation. The battery must be replaced shortly.
	Flashing: The battery must be replaced.
	Flashing in combination with **– – – –**: The battery is flat. The device switches itself off.

Operating SaneoTENS


GB

Prepare application

The device is switched off, the display blank.

➔ Connect the electrode cable with the connecting socket, channel A on the front side of the device.

➔ If you intend to use both channels:
Connect another electrode cable with the connecting socket channel B.

Prepare the skin areas in the application area:

Clean the area of skin to be treated from sweat, moisture and dirt. The skin must be dry (do not use any salves or creams) and undamaged.

➔ Take the body surface electrodes out of the protective packaging.

➔ Connect the electrodes with the electrodes cable/s (included in delivery).

Place the electrodes

When using 2 channels, if desired two different body parts can be treated. One electrode cable will be connected with two self-adhesive electrodes per channel.

➔ Place the electrodes according to the body region and the selected program (see chapter **Program overview by indication** on page 52 and the electrode placements shown on the **cover of this instruction manual**). Please observe the safety instructions (page 43 et seq.)!

➔ After treatment, the electrodes can be stuck back onto the backing foil and stored in the protective packaging.

Switching SaneoTENS on and off

SaneoTENS is switched on and off via the **On/Off** button. After the system has been switched on the software version will briefly be displayed. To switch off the **On/off** button must be pressed for approx. 3 seconds.

Set the program

To set the desired program (see chapter **Program overview**, page 50):

➔ Please switch on the SaneoTENS using the **on/off** button. The *Application menu* of the preset program will be shown briefly in the display.

➔ Hold the **Menu/confirmation** button down for about 3 seconds. The display now shows the *Program selection menu* with the current program.

➔ Press the ▲ or ▼ button until the desired program is displayed.

➔ Briefly press the **Menu/confirmation** button. The *Program selection menu* is completed. The selected program will immediately be displayed in the *Application menu*.

➔ Then select the individual stimulation intensity with the ▲/▼ buttons.

Please note: the program can only be changed when the device is not producing current. To change a program while it is running, the stimulation intensity must first be reduced to “0” using the ▼ button.

Select stimulation intensity

Monitor your own sensitivity!

The current feel or stimulation intensity must be clearly able to be felt, but in no way be uncomfortable or even painful.

The stimulation intensity for channel A or channel B will be displayed on the left or right in the *Treatment menu* (unit: **mA**).

Please select the stimulation intensity for **channel A** using the ▲/▼ buttons to the **left** of the display and for **channel B** using the ▲/▼ buttons to the **right** of the display.

➔ A short press of the ▲ button increases the stimulation intensity by 1 mA for the corresponding channel.

➔ A short press of the ▼ button reduces the stimulation intensity by 1 mA for the respective channel.

The treatment begins.

If a suitable stimulation intensity was set and not changed for 20 seconds, the **Lock** symbol appears in the display. The keypad is secured.

➔ If necessary, unlock the keypad with the ▼ button, e.g. to adjust the intensity.

Treatment frequency

As long as no side effects such as skin irritation occur, you can perform several treatments per day if necessary, ideally 3 times a day for 30 minutes each.

Programming freely available programs

The therapy programs **u1, u2** and **u3** can be set individually as required.

Program	Description
u1	Base program continuous stimulation (Frequency and program duration configurable)
u2	Base frequency modulation program (frequency, increase, plateau, decrease and pause time and program duration can be configured)
u3	Base program burst stimulation (pause time and program duration can be configured)

Parameter overview

Depending on the base program only specific parameters can be set. For the meaning of the individual parameters please see the following table. The value range of the parameter is given (in brackets).

Parameter	Description (value range)
Hz	Treatment frequency in Hz (2 to 120)
	Increase time in sec (0.5 to 5)
	Plateau time in sec (1 to 10)
	Decrease time in sec (0.5 to 5)
	Pause time in sec (1 to 10)
dEL	Pause between burst groups in sec (0.1 to 5)
	Treatment time in min (5 to 99)

Programming

The programming of medically effective, useful and combinable parameter values should only take place in consultation with your doctor or therapist.

Enter parameter

You can enter the settings of therapy programs **u1, u2** and **u3**. Your program configuration will then be called up like a preset program.

➔ Hold the **Menu/confirmation** button down. After about three seconds, the display shows the *program selection menu* with the current program.

➔ Press the ▲ or ▼ button until the therapy program you wish to configure is displayed.

➔ Press the **Menu/confirmation** button to navigate back to the first user definable parameters.

➔ The parameter value to be set flashes and can be decreased by pressing ▼ or increased by pressing▲.

➔ Press the **Menu/confirmation** button to confirm the selected value. The next parameter will be displayed.

➔ Continue to the last parameter **Treatment time** as described above.

➔ After you have also confirmed the treatment time with the **Menu/confirmation** button, the *Program selection menu* will end.

➔ The *Application menu* for the configured program appears in the display. The treatment with the previously set parameter value can begin.

Start the treatment

To start the treatment, set the stimulation intensity as described in the chapter **Select stimulation intensity** (page 62).

End treatment, cleaning

After the treatment time has finished, the device automatically switches itself off. The treatment can be interrupted at any time using the **On/off** button.

➔ Only remove the electrode(s) once the device is switched off.

➔ After use, carefully clean the electrodes and the device with a slightly damp, lint-free cloth. Cleaning the electrodes increases their service life.

➔ Stick the electrodes back on the carrier film and store the backing foil in the protective packaging.


If the adhesive force of the electrodes deteriorates, simply wet the adhesive surface with a bit of water and store the electrodes somewhere cool (e.g. in a refrigerator) until the next use. The electrodes must only be replaced when this no longer has any effect.

GB

Battery change

➔ Switch the device off before changing the battery.

1 Open the battery compartment.
2 Correctly insert the battery in the battery compartment.
3 Close the battery compartment.

➔ Please observe the applicable legal regulations in your country regarding the environmentally friendly disposal of batteries, accessories and packaging.


Technical data and support

Technical specifications

Name and model	SaneoTENS
Type	2-channel electrostimulation system
German Medical Devices Act (MPG) class	IIa
Dimensions	105 mm x 65 mm x 26 mm
Weight	Approx. 130 g including batteries
Power supply	Alkaline block battery or NiMH accumulator (9V)
Classification	BF
Current form	Biphasic rectangular pulse, charge-compensated
Operating conditions	Ambient temperature of +5 °C to +40 °C relative humidity of 15% to 90%, non-condensing, but without required water vapour partial pressure of more than 50hPa; air pressure from 700hPa to 1060hPa
Transport and storage conditions	Temperature range: -25°C to +5°C; relative humidity up to 90% at +5°C to 35°C (without condensation); > 35°C to 70°C at water vapour pressure up to 50hPa
IP classification:	IP22, protected against access with e.g. tools >=12.5mm and against dripping water
Operating time/number of applications in operation	20 applications (depending on the program)
Expected service life	15 years
Channels	2
Current strength	max. 60 mA ± 5 mA at 500 Ω
Max. Output voltage	100 V ±15 %
Treatment time	5 min. to 99 min. (adjustable in 1 min increments)
Impulse width	200 µs to 400 µs (adjustable in 1 µs increments)
Frequency range	2 Hz to 120 Hz (adjustable in 1 Hz increments)
Time required for cooling down or warming up from an unfavourable storage temperature to operating readiness at an ambient temperature of +20 °C:	
From -25 °C to +5 °C warming	2 hours
From +70 °C to +40 °C cooling	2 hours

GB

Labelling

Label packaging

Label device

Disposal / Maintenance

Environmentally friendly disposal

The national legal regulations applicable in your country for the environmentally friendly disposal of accessories, batteries and packaging must be observed.

This product must be disposed of separately from household waste in accordance with the applicable legal regulations. Hand it in at an official collection point. The raw material reserves are preserved due to separate collection and recycling and it is ensured that all regulations for the protection of health and the environment are observed when recycling the product.

Electrodes (Germany)

In Germany, all electrode types listed in the scope of delivery can currently be disposed of with normal household waste after use.

Maintenance

The system is designed and manufactured in such a way that maintenance is not required within the expected product life span (15 years). Each time the device is switched on, it performs self-tests to check that it is functioning correctly and can be used safely.

Fault and problem handling

SaneoTENS is a certified and approved medical device. Therefore, please report all serious incidents occurring in connection with the product to the manufacturer and the responsible State authority (in Germany the responsible authority is the Federal Institute for Drugs and Medical Devices (BfArM)).

If the suggested solutions do not solve the problem, please contact the manufacturer or your dealer.

1. **The device does not switch on when the On/off button is pressed.**
 - Make sure that the battery is correctly inserted and is making contact.
 - You may need to replace the battery.

2. **No stimulation can be felt.**
 - Press the On/off button to restart the device.
 - Check the placement of the electrodes. The electrodes may not overlap and must be fully stuck to the skin.
 - Increase the stimulation in stages
 - Check the cable connections between the devices and the electrodes.

- Try to improve the adhesive force of the electrodes by moistening. If the adhesive force of the electrodes is exhausted the current distribution is no longer guaranteed. Please replace the electrodes.

3. The current set jumps to 0 mA.

- The cable must be connected with the device and the electrodes and the electrodes must be placed correctly so that stimulation is possible.
- Try to restore the adhesive force of the electrodes by moistening. If the adhesive force of the electrodes is exhausted the current distribution is no longer guaranteed. Please replace the electrodes.

4. The battery symbol flashes and when the stimulation intensity is boosted the device switches off.

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